


Compound digital LCD biological microscope DMS-651



FEATURES

DMS-651 is for medical and health units, universities, research institutes to observe and research micro organisms, cells, bacteria, tissue culture, suspended substances, sediments and so on, especially for continuous observation of cells, bacteria in cultivation liquid. And the microscope can easily snapshot the growing and dividing process. It's widely used for cytology, parasites, oncology, immunology, genetic engineering, industrial microbiology, botany and so on.

LAS IMAGENES PUEDEN VARIAR DEL ORIGINAL; INFORMACIÓN SUJETA A CAMBIO SIN PREVIO AVISO

Cualquier duda o aclaración favor de llamarnos, estamos para SERVIRLE

México: [-52] 55-5300-4517, 55-53004271, 55-5312-2536

Querétaro: [-52] 442-340-0250, 442-340-0251, 442-193-5678

Puebla: [-52] 222-219-9999, 222-887-0114, 222-228-1633

www.cmsmetrology.com.mx

www.controlymedicion.com.mx

SPECIFICATION

Digital parts

Pixel	5Mega
Maximal Resolution	2560x1920
Video resolution	640x480
LCD screen	8 inches TFT screen, resolution is 800x600
Video output	AV output (NTSC/PAL)
Data output	USB 2.0
Menu language	Chinese/english
Storage	SD Card slot with optional memory cards (maximal storage capacity is 4G)
Scene mode	Standard/softness/vivid
Exposure mode	Auto exposure
Date Mode	year, month, day, hour, minute

Optical parts

Eye piece	High point, extra wide field eyepiece EW 10x/22	
Objective	LWD infinite plan objective	4x/0.1, WD 18mm
		40x/0.6, WD 2.6mm
	Infinite plan phase objective	PH10x/0.25, WD 10mm
		PH20x/0, WD 5.1mm
Nosepiece	Quintuple nosepiece	
Condenser	ELWD condenser NA 0.3 LWD 72mm (without condenser 150mm)	
Annulus	10x, 20x, 40x pahse annulus plate	
Stage	plain stage 160x250mm	
	Glass insert	
	Auxiliary stages 70x180mm	
Focusing	Coaxial coarse and fine adjustment, vertical objective movement, coarse stroke	
	37.7" per rotation, fine stroke 0.2mm per rotation	
Illumination	Halogen lamp 6V/30W	
Filter	Diameter 45mm, blue green and ground glass	

LAS IMAGENES PUEDEN VARIAR DEL ORIGINAL; INFORMACIÓN SUJETA A CAMBIO SIN PREVIO AVISO

Cualquier duda o aclaración favor de llamarnos, estamos para SERVIRLE

México: [-52] 55-5300-4517, 55-53004271, 55-5312-2536
 Querétaro: [-52] 442-340-0250, 442-340-0251, 442-193-5678
 Puebla: [-52] 222-219-9999, 222-887-0114, 222-228-1633