




Microscopio Metalurgico portatil



Cualquier duda o aclaración favor de llamarnos, estamos para SERVIRLE

México: [-52] 55-5300-4517, 55-53004271, 55-5312-2536

Querétaro: [-52] 442-340-0250, 442-340-0251, 442-193-5678

Puebla: [-52] 222-219-9999, 222-887-0114, 222-228-1633

Portable Trinocular Metallurgical Microscopes



Specification

Model No.	Microscope Magnification	X-Y Travel Range	Objective Lens
DP-PM9000	50X/100X/200X/500X	No Limit	5X/10X/20X/50X



Standard Module

- (1) 10X Trinocular Microscope Head with WF 10X/18 eyepiece
- (2) Coaxial Illumination Cast-iron Base
- (3) Plain achromatic convertible objective lens unit
- (4) CCD Connect Unit (CCD Optional)
- (5) X-Y-Z Cross moveable slide unit
- (6) 6V 20W Halogen light intensity adjustable illumination unit
- (7) Focusing coaxial coarse and fine adjustment unit
- (8) X: 30mm; Y: 49mm; Z: 70mm : 1mm
X axis double roller unit, moveablerange: no limit, fine adjustment range: 30mm; Y axis moveable range: 49mm; Z axis focus moveable range:70mm, fine adjustment range: 1mm.